

Eye Wash in 2-Pack CEDERROTH



REF 7251



Eye Wash



GB Eye Wash...
FR Double vision...
NL Oogwasmiddel...
IT Lavaggio Per...
ESPT Colirio...
CE 0413

REF 7251



Eye Wash



GB Eye Wash...
FR Double vision...
NL Oogwasmiddel...
IT Lavaggio Per...
ESPT Colirio...
CE 0413

Cederroth Eye Wash works to neutralise the eye's pH value, i.e. it has a neutralising effect on alkali or acid splashes, giving a quicker result than a normal sodium chloride solution. The effect is generally stronger on alkalis than acids. Shelf life of 4.5 years. If you have an accident with a substance that contains calcium (for example, cement, concrete or lime), you can safely rinse with Cederroth Eye Wash as it is 100% phosphate-free. This ensures that the eye is rinsed as effectively as possible.

Content:

2 (x 500 ml) sterile isotonic borate buffered sodium chloride solution without preservatives for single use. Bottle dimensions: Ø 6.5 x 23.5 cm. Shelf life 4.5 years

REF: 725200